

Self-releasing Roll Belt CHECKLIST



Facility		Contact	
Facility Name		Name	
Street Address		Title	
City/State/ZIP		Nursing Unit	
Instructor		Date	
Signature			

- 1. Open the belt and center it across the mattress with the waistband facing up. The hook fastening section will be located on the right side, and the long straps will hang off both sides of the bed.
- 2. Position the patient in bed with the belt at waist level.
- 3. Bring the belt around the patient's waist and secure with the hook and loop fastener. Connect the quick-release buckle. The belt must be snug, but not interfere with breathing. To check for proper fit, slide an open hand (flat) between the belt and the patient.
- 4. Attach the two long straps on each side to a movable part of the bed frame, at waist level, using quick-release ties or buckles.
- 5. Check that the straps are secure and will not change position, loosen, or tighten if the patient pulls on them, or if the bed is adjusted.

Rx ONLY

Before use, make sure to read the instructions accompanying the product. A copy of the instruction sheet is included with the product and can also be downloaded at www.posey.com.

Posey Company • 5635 Peck Road, Arcadia, CA91006-0020 USA • Phone: 1.800.447.6739 • Fax: 1.800.767.3933 • www.posey.com