

Omni Belt 4116PQ CHECKLIST



Facility		Contact	
Facility Name		Name	
Street Address		Title	
City/State/Zip		Nursing Unit	
Instructor		Date	
Signature			

Description	Competency Level	Date	Comments
Competency Level Definition - (1) Needs remediation and re-evaluation, (2) Performs skill with accuracy under supervision, (3) Performs skill with accuracy independently			
APPLICATION INSTRUCTIONS: FOR WHEELCHAIR OR GERI-CHAIR:			
1. Open the lap belt and position it so that the patient will sit on it.			
2. Bring the ends of the connecting straps down at a 45-degree angle between the seat and wheelchair sides or around the outside of the chair frame (Geri-chair).			
3. Criss-cross the straps behind the chair and draw them around the opposite side kick spurs (wheelchair), or the lower cross braces (Geri-chair). Insert the strap ends through the white "D" ring on each side, and attach the hook and loop fasteners.			
4. Use slide buckles to remove any slack. Check that the straps are secure and will not change position, loosen, or tighten if the patient pulls on them, or if the chair is tilted or adjusted.			
5. Position the patient as far back in the seat as possible, with the buttocks against the back of the chair.			
6. Lay the lap belt across the patient's thighs and connect the quick release buckle so the belt is across the patient's lower lap.			

Rx ONLY

Before use, make sure to read the instructions accompanying the product. A copy of the instruction sheet is included with the product and can also be downloaded at www.posey.com.

5635 Peck Road, Arcadia, CA 91006-0020 USA • Phone: 1.800.447.6739 • Fax: 1.800.767.3933 • www.posey.com



Description	Competency Level	Date	Comments
Competency Level Definition - (1) Needs remediation and re-evaluation, (2) Performs skill with accuracy under supervision, (3) Performs skill with accuracy independently			
7. Use quick-release buckle to adjust the lap belt. The belt must be snug, but not interfere with breathing. To check for proper fit, slide an open hand (flat) between the belt and the patient.			
FOR HOSPITAL BED USE: 1. Open the belt and center it across the mattress.			
2. Position the patient in bed with the belt at waist level.			
3. Bring the belt around the patient's waist and connect the quick release buckle.			
4. Use quick-release buckle to adjust the belt. The belt must be snug, but not interfere with breathing. To check for proper fit, slide an open hand (flat) between the belt and the patient.			
5. Feed the end of each connecting strap around a movable part of the bed frame, at waist level. Insert the strap ends through the white "D" ring on each side, and attach the hook and loop fasteners.			
6. Use slide buckles to remove any slack. Check that the straps are secure and will not change position, loosen, or tighten if the patient pulls on them, or if the bed is adjusted.			