

Quick-Release Twice-As-Tough™ Cuffs 2789Q CHECKLIST



Facility		Contact	
Facility Name		Name	
Street Address		Title	
City/State/Zip		Nursing Unit	
Instructor		Date	
Signature			

Description	Competency Level	Date	Comments
Competency Level Definition - (1) Needs remediation and re-evaluation, (2) Performs skill with accuracy under supervision, (3) Performs skill with accuracy independently			
2789Q Wrist, pair, single strap with quick-release buckles (color-coded blue)			
<p>APPLICATION INSTRUCTIONS:</p> <p>⚠ WARNING</p> <p>Connecting Straps Before each use, check cuffs and straps for cracks, tears, and/or excessive wear or stretch; cracked or broken buckles or locks; and/or that hook and loop adheres securely, as these may allow patient to remove cuff. Discard if device is damaged.</p> <p>If you have any questions about patient safety, ask the appropriate medical authority for alternatives.</p>			
1. Follow these steps to apply device: a. Attach the female end of the quick-release buckle (short strap) to a movable part of the bed frame, out of the patient's reach. Secure by feeding the female end through the loop in the strap.			
b. Insert the male end of the connecting strap into the female end of the short strap. Listen for a "snapping" sound. Pull firmly on straps to ensure a good connection.			

Rx ONLY

Before use, make sure to read the instructions accompanying the product. A copy of the instruction sheet is included with the product and can also be downloaded at www.posey.com.

Posey Company • 5635 Peck Road, Arcadia, CA 91006-0020 USA • Phone: 1.800.447.6739 • Fax: 1.800.767.3933 • www.posey.com



Description	Competency Level	Date	Comments
Competency Level Definition - (1) Needs remediation and re-evaluation, (2) Performs skill with accuracy under supervision, (3) Performs skill with accuracy independently			
Applying the Cuffs (repeat steps 1-4 for each limb) 1. Wrap the neoprene piece (the blue side should be positioned against the skin) around the wrist. Attach the black hook and loop pieces together, followed by the blue hook and loop pieces. The fuzzy piece should be sandwiched between the two pieces of hook. Be sure to overlap at least one inch (3 cm).			
2. Press the hook and loop closure together firmly and make sure it adheres securely. Slide ONE finger (flat) between the cuff and the inside of the patient's wrist/ankle to ensure proper fit. The cuffs must be snug enough to prevent escape, but not interfere with circulation.			
3. Release the quick-release buckle, twist, and reconnect. Listen for a "snapping" sound. Pull firmly on straps to ensure a good connection.			
4. Adjust the bed strap(s) to allow desired freedom of movement, without compromising patient or caregiver safety.			