



Posey Company

5635 Peck Road

Arcadia, CA 91006-0020 USA

Toll Free: (800) 44-POSEY or (800) 447-6739

Web: www.posey.com

FAQ No.: 10

FAQ:

How many times can a Posey product be washed and maintain its product integrity?

Posey products should always be washed and dried according to the instructions that accompany the product. The life of the product depends on the number of times washed, wash cycle duration, detergent usage, bleach usage, wash water temperature, and the drying cycle temperature and duration.

Posey products should always be inspected before use for broken stitches or parts, torn, cut or frayed straps or fabric, or hook and loop fasteners that do not hold securely. If these components fail, it could result in injury or death. Destroy and dispose of products worn beyond their useful life.

For complete cleaning and use instructions for all Posey products, please refer to package insert instructions accompanying each product. Copies of all package insert instructions are available on the Posey website at www.posey.com.

We hope that this answers your questions. If you have any further questions, please feel free to call Posey Customer Service at (800) 44-Posey or (800) 447-6739.