



Posey Company

5635 Peck Road

Arcadia, CA 91006-0020 USA

Toll Free: (800) 44-POSEY or (800) 447-6739

Web: www.posey.com

FAQ No.: FAQ 17

FAQ: What are the weight limitations of Orthotic Drop Seat (#6331) and Drop Seat (#6540)?

The Posey Drop Seat and Orthotic Drop Seat are designed to withhold the weight capacity of the wheelchair for which they are intended. Although most standard wheelchairs are designed to withstand from 250 to 350 pounds, please check with the wheelchair manufacturer for the specific weight limit of the brand of chair you are using. Always install an anti-tipping device when using a Posey Drop Seat.

For complete cleaning and use instructions for all Posey products, please refer to package insert instructions accompanying each product. Copies of all package insert instructions are available on the Posey website at www.posey.com.

We hope that this answers your questions. If you have any further questions, please feel free to call Posey Customer Service at (800) 44-Posey or (800) 447-6739.