



Posey Company

5635 Peck Road

Arcadia, CA 91006-0020 USA

Toll Free: (800) 44-POSEY or (800) 447-6739

Web: www.posey.com

FAQ No.: 36

FAQ:

What are the Warranty and Cleaning Instructions for 8283 Posey Disposable Bed Sensor?

What is the length of warranty on Product 8283 Over the Mattress Alarm Pad?

- The warranty is 30 days as written on the product label.

Is the product Single Patient Use or can it be used from Patient to Patient ?

- The product labels states "Single Patient Use."

What are the cleaning instructions ?

- On the product label, it states to wipe-clean with liquid disinfectant.

Note: Do not submerge in liquid or sterilize with heat.

For complete cleaning and use instructions for all Posey products, please refer to package insert accompanying each product.

We hope that this answers your questions. If you have any further questions, please feel free to call Posey Customer Service at (800) 44-Posey or (800) 447-6739.