



# Posey Company

5635 Peck Road

Arcadia, CA 91006-0020 USA

Toll Free: (800) 44-POSEY or (800) 447-6739

Web: [www.posey.com](http://www.posey.com)

**FAQ No.:** 42

**FAQ:** Can Posey Limb Restraints be applied with one arm positioned upward, and the other arm positioned downward ?

---

Posey does not recommend the application of limb holder to restrain one arm over the head and the other arm restrained by the patient's side.

Complications and injury could occur. These complications could include, but are not limited to undue pressure in the axilla (armpit) region, and possibly shoulder dislocation to the arm that is positioned above the head, if the patient is aggressive.

For Safety, Application, and Cleaning instructions for all Posey products, please refer to package insert instructions accompanying each product. Copies of all package insert instructions are available on the Posey website at [www.posey.com](http://www.posey.com).

We hope that this answers your questions. If you have any further questions, please feel free to call Posey Customer Service at (800) 44-Posey or (800) 447-6739.