



Posey Company

5635 Peck Road

Arcadia, CA 91006-0020 USA

Toll Free: (800) 44-POSEY or (800) 447-6739

Web: www.posey.com

FAQ No.: 44

FAQ: Can Posey Mitts and Limb Holders be used in Wheelchairs & Chairs ?

Yes, under the following conditions:

Do not use Posey products on Toilets or any type of furniture, which does not allow application as directed in the product application sheet.



Make sure straps are secured at a juncture of the frame and will not slide in any direction, changing position of the device.

Make sure the patient is not able to slide down, or fall off a chair seat, possibly causing a dislocated shoulder.

Please refer to "Selecting The Right Posey Product" Page, in the Posey Healthcare Product Selection Guide, for products to reduce sliding in chairs.

For Safety, Application, and Cleaning instructions for all Posey products, please refer to package insert instructions accompanying each product. Copies of all package insert instructions are available on the Posey website at www.posey.com.

We hope that this answers your questions. If you have any further questions, please feel free to call Posey Customer Service at (800) 44-Posey or (800) 447-6739.