



Posey Company

5635 Peck Road

Arcadia, CA 91006-0020 USA

Toll Free: (800) 44-POSEY or (800) 447-6739

<http://www.posey.com>

FAQ No.: 47

FAQ: The Posey 8070 bed uses a 2-prong plug. Does this plug safely protect the user from electrical shock like a grounded 3-prong plug?

The Posey 8070 bed uses a method called “double insulation” rather than “grounding” to protect the user from electrical shock. Both grounded and double insulated systems are designed to reduce the risk of electrical shock to the user. The grounded system accomplishes this by providing a ground path for any electrical fault that may occur in the appliance, thus preventing the user from becoming the ground path. The double-insulated system, because it contains no path to ground and isolates the user from any electrical faults inside the appliance, accomplishes this in a different manner.

Both systems are recognized and approved by all national listing bodies, *i.e.*, UL, CSA, TUV, as well as by the United States NFPA, NEC, and OSHA. The Posey 8070 bed is a UL Listed Medical Equipment with UL File Number E37927.

We hope that this answers your questions. If you have any further questions, please feel free to call us at 800-44-POSEY.

For complete cleaning and use instructions for all Posey products, please refer to package insert instructions accompanying each product. Copies of all package insert instructions are available on the Posey website at <http://www.posey.com>.

We hope that this answers your questions. If you have any further questions, please feel free to call Posey Customer Service at (800) 44-Posey or (800) 447-6739.